

BloodStor[®] 100

FREEZE MEDIA

BEST-IN-CLASS, OPTIMIZED BIOPRESERVATION MEDIA FOR CELLS AND TISSUES

100% DMSO

Sterile, USP Vial

**Standard and Custom
Fill Volumes**

**USP/Highest Quality
Components**

Filled Under cGMP

**Sterility and
Endotoxin Testing**



BLOODSTOR[®]

BloodStor[®] 100 is a series of cGMP filled DMSO freeze media products used to cryopreserve stem cells and other cells isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor products are formulated using USP grade components. Standard and custom packaging and fill volumes are available.

ORDERING INFORMATION

Product Name	Size	Part #	Retail Price
BloodStor [®] 100	1 x 50mL vial	410301	\$32
BloodStor [®] 100	1 x 100mL vial	410302	\$60

RELATED PRODUCTS

Product Name	Size	Part #	Retail Price
BloodStor [®] 55-5	16 x 5.6mL vials per case	355371	\$448
BloodStor [®] 55-5	16 x 7.2mL vials per case	355372	\$528

MATERIALS ARE MANUFACTURED UNDER cGMP

TEST	METHOD	LIMITS
Visual Inspection	Visual Inspection	Clear to slightly yellow solution with no visible particulates
pH	SOP 3006	Report Results
Endotoxin	Kinetic Chromogenic USP <85>	≤ 1 EU/mL
Sterility	Membrane Filtration USP <71>	Sterile